



PYROIS
TECH

COB-series User Manual

Fiber coupled
LED light source
UV-VIS-NIR



Pyroistech S.L.

Campus Arrosadia, Ed. Los Tejos,
Pamplona 31006, Navarra, Spain

© 2019 Pyroistech

www.pyroistech.com
info@pyroistech.com

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, by any means, electronic, mechanical, photocopying, recording, or otherwise, without written permission from Pyroistech S.L.

This manual accompanies a product order and is subjected to the condition that it shall not, by way of trade or otherwise, be lent, re-sold, hired out or otherwise circulated without the prior consent of Pyroistech S.L. in any form of binding or cover other than that in which it is published.

Disclaimer

Every effort has been made to make this manual as complete and as accurate as possible, but no warranty is implied on it. The information provided is on an "as is" basis. Pyroistech S.L. shall have neither liability nor responsibility to any person or entity with respect to any loss or damages arising from the information contained in this manual.

Table of Contents

0. Introduction.....	5
1. Technical Specifications.....	5
2. Operation.....	6
2.1 Package Content.....	6
2.2 Before use.....	6
2.3 Maintenance.....	6
2.4 COB-series Diagrams.....	7
2.5 Operation Modes.....	8
3. Important Safety Notes.....	9
4. Warranty.....	10
5. Compliance.....	10

0. Introduction

Thank you for purchasing this COB light source from Pyroistech.SL. This document describes the COB light source and provides you with instructions for its correct operation. Please do not hesitate to contact us through info@pyroistech.com if you have any questions or doubts about this manual.

The COB series comprises high power LED light sources with peak wavelengths that range from UV to NIR (270 – 940 nm). They are conceived as user friendly, compact and silent tools, being ideal for fluorescence, spectroscopy and general optical fiber illumination applications.

These light sources have been designed to obtain a high coupling efficiency to optical fiber through a selectable SMA/FC connector. In addition, they can be operated either in continuous (CW) or pulse width modulation (PWM) modes.

1. Technical Specifications

Optical characteristics

Reference	λ (nm)	Power ¹ (mW)	FWHM (nm)	Reference	λ (nm)	Power ¹ (mW)	FWHM (nm)
COB-270	270	0.12	15	COB-460	460	10.7	24
COB-280	280	0.11	15	COB-523	523	4.8	36
COB-310	310	0.50	15	COB-590	590	2.0	15
COB-365	365	11.8	12	COB-623	623	10.3	17
COB-385	385	10.6	12.5	COB-660	660	10.6	18
COB-395	395	11.0	16	COB-740	740	7.3	18
COB-405	405	11.4	20	COB-840	840	13.1	33
COB-457	457	13.2	20	COB-940	940	29.0	40
COB-6500 White	-	16.5	-	COB-3000 White	-	2.4	-

(1) Typical output power measured with an optical fiber (core diameter 600 μ m) and a Si photodiode (OPHIR™)

Electrical characteristics

Power input	V_{IN}	10 - 18 V
	I_{IN}	1.5 A max
	Connector type	DC female 2.1mm
Signal input for pulsed mode	V_{IN}	5V
	I_{IN}	5 mA typ.
	Connector type	SMA female

Other characteristics

Output connector	SMA 905 or FC
Working T	0 - 50 °C
Humidity	< 80% RH
Equipment Surface T	40°C
Stabilization time	25 min typ.
Size	10.5 x 11.5 x 8 cm
Weight	500 g

2. Operation

2.1 Package Content

- Requested COB light source
- Power source: input 100 - 240 V~, 50 - 60 Hz; output 18 V, 2 A; 99 x50 x33 cm, cord length 1.5 m.
- Power cord 1.85 m long
- User's manual

Inspect carefully the devices and make sure there is no damage. On the contrary, do not employ the light source and please contact us through info@pyroistech.com for repair or replacement information.

2.2 Before use

Insert the plug of the power cord into the power supply and connect the power cord to the power outlet. Then, introduce the plug of the power supply into the light source power input (DC female type).

Unscrew the cap of the SMA 905/FC connector on the front face of the light source and connect the fiber cable. Always do this step before turning on the light source. Now, the COB light source is ready to be used.

The COB source is designed to have intensity control in continuous and modulation modes. The easy user interface permits to select the desired mode (see section [2.4](#)).

2.3 Maintenance

The output connector of the light source should be periodically cleaned employing compressed air to remove the dust. It is recommended to have the output connector with the cap on it whenever the light source is not being used. The connector of the fiber used with the light source also has to be clean to ensure the best possible performance.

2.4 COB-series diagrams

1	Intensity regulator (20% - 100%) in CW mode
2	LED indicator
3	Switch (ON-OFF-PWM)
4	PWM input (5V max)
5	Power input (10 – 18 V, 1.5 A max)
6	SMA 905/FC output connector



2.5 Operation Modes

The controls feature a three way switch (3) to select the operation mode: continuous, PWM or off. A second control (1) is used to manually adjust the drive current in continuous mode. A connector allows (4) to introduce a control signal in the PWM mode.

• OFF mode

The switch (3) on the right of the frontal face of the light source has to be on the middle position ('OFF'). It is suggested to have the switch in this position when connecting the power supply to the light source through (5).

• CW mode

To select the continuous mode, please move the switch (3) to the left position ('ON'). In this mode, the left knob (1) controls the intensity. The intensity can be adjusted between 20% (fully counter-clockwise) and 100% (fully clockwise) of maximum intensity.

• PULSED mode

To work in pulsed mode, please move the switch (3) to the right position (pulsed signal drawing). In this mode, the intensity can be adjusted between 0% and 100% of maximum intensity changing the duty cycle of the PWM signal introduced through the PWM input connector (4) of SMA female type.

The PWM signal must vary between 0 and 5 V. Voltages out of this range can damage the PWM input and disable the PWM mode. Voltages below 1.5 V will be recognized as low level while tensions above 2.8 V will be considered high level. Pulsed frequencies up to 100 kHz can be used.

Whenever the LED is on, either in CW mode or in pulsed mode; the green LED indicator (2) located on the front face of the light source will be on. In case of being in pulsed mode, the LED indicator will blink in the same way as the light source.

The light comes out through a connector (6) that can be SMA 905 or FC depending on the customer's choice and that has been specially designed to ensure that the optical fiber can be correctly screwed to the light source.

3. Important Safety Notes

- Do not remove or alter any installed safety device on this equipment. Doing so will cancel your warranty and create an unsafe operating environment.
- There are NO user serviceable parts inside. Dangerous currents are present in this device. Only allow qualified personnel to service this unit.
- Inspect this unit and its power supply before using it for the first time. Do not use the unit if it is damaged in any way. Contact us for repair or replacement information.
- During operation do not cover the source or obstruct the air flow for its refrigeration. Avoid exposure to direct sun light. A rise in the light source's temperature could affect its operation or even damage its components.
- The output connector of the light source may get hot during operation. After its employment, allow enough time to cool down before handling.
- Optical radiation can damage your eyes. Do NOT stare directly at the light beam.



CAUTION !

Do NOT stare directly at the light beam

- Proper protective eyewear must be worn when using light sources that emit UV radiation (COB-270, COB-280, COB-310, COB-365, COB-385, COB-395, COB-405). Avoid exposure to the beam. It is hazardous to skin and eyes ,and may cause cancer.



WARNING ! UV LEDs

Avoid eye and skin exposure to the emitted UV light

- COB-840 and COB-940 emit non visible infrared light, which can be hazardous depending on total system configuration (including, but not limited to optics, drive current and temperature). Observe safety precaution given in IEC 62471 when operating these light sources.

4. Warranty

The COB light source is covered by Pyroistech's 2 years warranty.

5. Compliance

This device has been tested and complies with the following standards:

EMC IEC 61000-4-2:2008 IEC 61000-4-8:2009
IEC 61000-4-3:2006 CISPR 11:2009
IEC 61000-4-4:2012 Directive 2014/30/EU
IEC 61000-4-5:2014 EN 61326-1:2013
IEC 61000-4-6:2013

Safety IEC 61010-1:2010



EMC 2014/30/EU
RoHS-compliant



WEEE
Compliance

Version	Date	Description
1.0	9 September 2019	First Document

PYROISTECH S.L.

Campus Arrosadía, Los Tejos, IEEC

31006, Pamplona, Spain

+34 948 168 845

Email: info@pyroistech.com



PYROIS
TECH

"riding the light"

www.pyroistech.com



<https://es.linkedin.com/company/pyroistech>