



PYROIS
TECH

COB-series User Manual

Fiber coupled
LED light source
UV-VIS-NIR



Pyroistech S.L.
C/ Tajonar 22
Edificio Jerónimo de Ayanz
31006, Pamplona, Spain
+34 698 91 04 12

www.pyroistech.com
info@pyroistech.com

© 2023 Pyroistech

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, by any means, electronic, mechanical, photocopying, recording, or otherwise, without written permission from Pyroistech S.L.

This manual accompanies a product order and is subjected to the condition that it shall not, by way of trade or otherwise, be lent, re-sold, hired out or otherwise circulated without the prior consent of Pyroistech S.L. in any form of binding or cover other than that in which it is published.

The specifications indicated in this manual are subject to change without prior notice.

Disclaimer

Every effort has been made to make this manual as complete and as accurate as possible, but no warranty is implied on it. The information provided is on an "as is" basis. Pyroistech S.L. shall have neither liability nor responsibility to any person or entity with respect to any loss or damages arising from the information contained in this manual.

Table of Contents

0. Introduction.....	5
1. Technical Specifications.....	5
2. Operation.....	6
2.1 Package Content.....	6
2.2 Before use.....	6
2.3 Maintenance.....	6
2.4 COB-series Diagrams.....	7
2.5 Operation Modes.....	8
3. Important Safety Notes.....	9
4. Warranty.....	10
5. Compliance.....	10

0. Introduction

Thank you for purchasing this COB light source from Pyroistech.SL. This document describes the COB light source and provides you with instructions for its correct operation. Please do not hesitate to contact us through **info@pyroistech.com** if you have any questions or doubts about this manual.

The COB series comprises high power LED light sources with peak wavelengths that range from UV to NIR (270 – 1050 nm). They are conceived as user friendly, compact and silent tools, being ideal for fluorescence, spectroscopy and general optical fiber illumination applications.

These light sources have been designed to obtain a high coupling efficiency to optical fiber through a selectable SMA/FC connector. In addition, they can be operated either in continuous (CW) or pulsed modes.

1. Technical Specifications

Fiber coupled LED Light Source- Ultraviolet (UV)

Optical characteristics			
Name	Peak wavelength	Spectral range FWHM	Typical optical power output (2)
COB-265	265 nm	11 nm	TBD
COB-270	270 nm	15 nm	0.12 mW
COB-280	280 nm	15 nm	0.11 mW
COB-300	300 nm	20 nm	TBD
COB-310	310 nm	15 nm	0.50 mW
COB-325	325 nm	12 nm	TBD
COB-340	340 nm	10 nm	TBD
COB-365	365 nm	12 nm	11.8 mW
COB-385	385 nm	12.5 nm	10.6 mW
COB-395	395 nm	16 nm	11.0 mW
COB-405	405 nm	20 nm	11.4 mW

(*) Measured with an optical fiber with a core diameter of 600 μ m and 0.22 NA

Fiber coupled LED Light Source- Visible (VIS)

Optical characteristics			
Name	Peak wavelength	Spectral range FWHM	Typical optical power output (2)
COB-430	430 nm	17 nm	TBD
COB-457	457 nm	20 nm	13.2 mW
COB-460	460 nm	24 nm	10.7 mW
COB-523	523 nm	36 nm	4.8 mW
COB-590	590 nm	15 nm	2.0 mW
COB-623	623 nm	17 nm	10.3 mW
COB-660	660 nm	18 nm	10.6 mW
COB-EX White	EX White	TBD	0.75 mW
COB-3000 White	3000 White	TBD	2.4 mW
COB-6500 White	6500 White	TBD	16.5 mW

(*) Measured with an optical fiber with a core diameter of 600 μ m and 0.22 NA

Fiber coupled LED Light Source- Infrared (IR)

Optical characteristics			
Name	Peak wavelength	Spectral range FWHM	Typical optical power output (2) (*)
COB-740	740 nm	18 nm	7.3 mW
COB-840	840 nm	33 nm	13.1 mW
COB-940	940 nm	40 nm	29.0 mW
COB-1050	1050 nm	55 nm	58.1 mW (**)
COB-1100	1100 nm	50 nm	TBD
COB-1200	1200 nm	65 nm	TBD
COB-1300	1300 nm	80 nm	TBD
COB-1450	1450 nm	95 nm	TBD
COB-1550	1550 nm	102 nm	TBD
COB-1650	1650 nm	120 nm	TBD

(*) Measured with an optical fiber with a core diameter of 600 μ m and 0.22 NA

(**) Measured with a VIS-NIR optical fiber with a core diameter of 1000 μ m and 0.5 NA.

Electrical characteristics

Power input	V_{IN}	10 - 18 V
	I_{IN}	1.5 A max
	Connector type	DC female 2.1mm
Signal input for pulsed mode	V_{IN}	5V
	I_{IN}	5 mA typ.
	Connector type	SMA female

Other characteristics

Output connector	SMA 905 or FC
Working T	0 - 50 °C
Humidity	< 80% RH
Equipment Surface T	40 °C
Stabilization time	25 min typ.
Size	10.5 x 11.5 x 8 cm
Weight	500 g

2. Operation

2.1 Package Content

- Requested COB light source
- Power source: input 100 - 240 V~, 50 - 60 Hz; output 18 V, 2 A; 99 x 50 x 33 cm, cord length 1.5 m.
- Power cord 1.85 m long
- User's manual

Inspect carefully the devices and make sure there is no damage. On the contrary, do not employ the light source and please contact us through info@pyroistech.com for repair or replacement information.

2.2 Before use

Insert the plug of the power cord into the power supply and connect the power cord to the power outlet. Then, introduce the plug of the power supply into the light source power input (DC female type).

Unscrew the cap of the SMA 905/FC connector on the front face of the light source and connect the fiber cable. Always do this step before turning on the light source. Now, the COB light source is ready to be used.

The COB source is designed to have intensity control in continuous and modulation modes. The easy user interface permits to select the desired mode (see section **2.4**).

2.3 Maintenance

The output connector of the light source should be periodically cleaned employing compressed air to remove the dust. It is recommended to have the output connector with the cap on it whenever the light source is not being used. The connector of the fiber used with the light source also has to be clean to ensure the best possible performance.

2.4 COB-series diagrams

1	Intensity regulator (20% - 100%) in CW mode
2	LED indicator
3	Switch (ON-OFF-PWM)
4	PWM input (5V max)
5	Power input (10 – 18 V, 1.5 A max)
6	SMA 905/FC output connector



2.5 Operation Modes

The controls feature a three way switch (3) to select the operation mode: continuous, pulsed or off. A second control (1) is used to manually adjust the drive current in continuous mode. A connector allows (4) to introduce a control signal in the PWM mode.

• OFF mode

The switch (3) on the right of the frontal face of the light source has to be on the middle position ('OFF'). It is suggested to have the switch in this position when connecting the power supply to the light source through (5).

• CW mode

To select the continuous mode, please move the switch (3) to the left position ('ON'). In this mode, the left knob (1) controls the intensity. The intensity can be adjusted between 20% (fully counter-clockwise) and 100% (fully clockwise) of maximum intensity.

• PULSED mode

To work in pulsed mode, please move the switch (3) to the right position (pulsed signal drawing). In this mode, the intensity can be adjusted between 0% and 100% of maximum intensity changing the duty cycle of the PWM signal introduced through the PWM input connector (4) of SMA female type.

The PWM signal must vary between 0 and 5 V. Voltages out of this range can damage the PWM input and disable the PWM mode. Voltages below 1.5 V will be recognized as low level while tensions above 2.8 V will be considered high level. A pulsed frequency of 500 Hz or lower is recommended to achieve the maximum resolution.

Whenever the LED is on, either in CW mode or in pulsed mode; the green LED indicator (2) located on the front face of the light source will be on. In case of being in pulsed mode, the LED indicator will blink in the same way as the light source.

The light comes out through a connector (6) that can be SMA 905 or FC depending on the customer's choice and that has been specially designed to ensure that the optical fiber can be correctly screwed to the light source.

3. Important Notes

Before operation

- Do not remove or alter any installed safety device on this equipment. Doing so will cancel your warranty and create an unsafe operating environment.
- There are NO user serviceable parts inside. Dangerous currents are present in this device. Only allow qualified personnel to service this unit.
- Inspect this unit and its power supply before using it for the first time. Do not use the unit if it is damaged in any way. Contact us for repair or replacement information.

During operation

- Do not cover the source or obstruct the air flow for its refrigeration. Avoid exposure to direct sun light. A rise in the light source's temperature could affect its operation or even damage its components.
- The output connector of the light source may get hot during operation. After its employment, allow enough time to cool down before handling.
- The equipment should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment should be observed to verify normal operation in the configuration in which it will be used.

Eye safety

- Optical radiation can damage your eyes. Do NOT stare directly at the light beam.



CAUTION !

Do NOT stare directly at the light beam

- Proper protective eyewear must be worn when using light sources that emit UV radiation (COB-270, COB-280, COB-310, COB-365, COB-385, COB-395, COB-405). Avoid exposure to the beam. It is hazardous to skin and eyes, and may cause cancer.



WARNING ! UV LEDs

Avoid eye and skin exposure to the emitted UV light

- COB-840 and COB-940 emit non visible infrared light, which can be hazardous depending on total system configuration (including, but not limited to optics, drive current and temperature). Observe safety precaution given in IEC 62471 when operating these light sources.

Electromagnetic Compatibility

- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of this equipment. Otherwise, degradation of the performance of this equipment could result.
- The use of cables, power supplies, transformers and accessories other than those specified by Pyroistech S.L., in addition to not being covered by the warranty, may lead to increased emissions and/or decreased electromagnetic compatibility.

4. Warranty

Pyroistech's general warranty for a product has a 1 year duration.

This warranty includes repairs and replacement of damaged parts due to a malfunction of the source, as long as said malfunction can be attributed to errors made in the manufacture of it by Pyroistech SL, not to a punctual misuse of the source or to a continued incorrect employment of it by the user, whether conscious or unconscious, due to not having followed the operation recommendations indicated by Pyroistech SL

There is the possibility of extending this guarantee. For more information, contact Pyroistech S.L.

5. Compliance

This device complies the following standards:



EMC 2014/30/EU
RoHS-compliant



Federal
Communications
Commision

Contact Pyroistech S.L. if you require more information about the electromagnetic compatibility of the product.



WEEE
Compliance

If you consider that the product has reached the end of its useful life and you want to dispose of it, you can contact Pyroistech S.L. so that it is in charge of its management.

Version	Date	Description
1.0	9 September 2020	First Document

PYROISTECH S.L.

C/ Tajonar 22

Edificio Jerónimo de Ayanz

31006, Pamplona, Spain

+34 698 91 04 12

Email: info@pyroistech.com



PYROIS
TECH

"riding the light"

www.pyroistech.com



<https://es.linkedin.com/company/pyroistech>